
Basic Cyanobacteria ID Guide



Jen Maucher Fuquay Phytoplankton Monitoring Network



DEFINITIONS

Cyanobacteria anatomy


Aerotope- gas vacuoles found in some cyanobacteria. Allow for change in position in water column based on light level and photosynthesis. Appear as clear refractive clusters with irregular outlines in vegetative cells. Best viewed with bright field microscopy.

Akinete- specialized cells for survival of adverse environmental conditions; resting

stage. Can germinate into viable vegetative cells when conditions improve. Akinetes are thick walled and appear granular.

Heterocyte (heterocyst)- site of nitrogen fixation. Non-vegetative cells that are transparent and have thick walls which may give 'halo' like appearance.

Vegetative cell- somatic (non-reproductive) cell




Filament- trichome and sheath

Sheath– a mucilaginous outer layer that encapsulates entire trichome


Intercalary- in between other cells i.e. not at terminal end of trichome

Facultative – optional or discretionary e.g. some cyanobacteria may or may not have aerotopes at any given time in their life cycle.



Vegetative cell


NO SO ES ORA



Akinete

APHANIZOMENON

General Characteristics




Straight, unbranched trichomes Trichomes tend to be in parallel bundles (AKA 'grass clippings') No sheath Order Nostocales

APHANIZOMENON




Things to look for:

- A. Aerotopes (facultative)
- B. Heterocytes are intercalary (facultative)
- C. Akinetes usually cylindrical and intercalary
- D. Terminal ends are elongated and may be "empty" looking



DOLICHOSPERMUM

(formerly Anabaena)




Unbranched trichomes Trichomes can be straight, coiled or curvy End cells look similar to other vegetative cells Order Nostocales

DOLICHOSPERMUM



Things to look for:

- A. Cells rounded or barrel shaped with aerotopes
- B. Heterocytes are intercalary
- C. Akinetes are intercalary
- D. Cells constricted at cross walls



RAPHIDIOPSIS

General Characteristics




Unbranched trichomes Trichomes can be straight, coiled or curvy Heterocytes ALWAYS terminal Order Nostocales


RAPHIDIOPSIS



Things to look for:


- A. Cells cylindrical with aerotopes (facultative)
- B. Heterocytes (when present) are always terminal at one or both ends
- C. Akinetes (when present) usually 1-3 cells back from heterocytes
- -Terminal cells attenuate (get narrower)
- -Difficult to ID when lacking heterocytes





PLANKTOTHRIX

Formerly Oscillatoria




Unbranched trichomes Trichomes can be straight or wavy NO heterocytes NO akinetes Order Oscillatoriales

PLANKTOTHRIX

Things to look for:


A.Cells mostly shorter than wide or isodiametric B. LOTS of aerotopes throughout cells Terminal cells can be: C. narrowly conical D. broadly rounded or swollen (calyptrate)- not shown

- <u>No</u> heterocytes (not a N₂ fixer)
- <u>No</u> akinetes



MICROCYSTIS

General Characteristics



Colonial Colonies surrounded by mucilage Shape of colony & mucilage can help to ID NO heterocytes NO akinetes Order Choococcales

MICROCYSTIS



Things to look for:

A. Rounded cells with aerotopes

Cells in colony may be

- A. loosely associated
- B. clathrate
- C. densely packed

Mucilage can vary in thickness

